

# Anti-FE65 Rabbit antibody



Catalog Number: 168029

## Product name

Anti-FE65 Rabbit antibody

## Specificity

Human

## Antibody description

FE65 Rabbit monoclonal antibody

## Preparation

Antigen: A synthetic peptide of human FE65

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,  
40%Glycerol, 0.01% sodium azide and 0.05% BSA

## Storage

Store at 4°C short term. Store at -20°C long term.  
Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications


WB, IHC-P, ICC/IF

## Dilutions

WB: 1/1000

IHC: 1/201/20

## Validations



Western blot detection of FE65 in A549,U251 cell lysates using FE65 Rabbit mAb(1:1000 diluted).Predicted band size:77kDa.Observed band size:100kDa.