

Anti-Phospho-TrkA (tyr490)/TrkB (Tyr516) Rabbit antibody



Catalog Number: 169029

Product name

Monoclonal

Anti-Phospho-TrkA (tyr490)/TrkB (Tyr516) Rabbit antibody

Ig Type

Rabbit IgG

Specificity

Human

Applications

WB

Antibody description

Phospho-TrkA (tyr490)/TrkB (Tyr516) Rabbit monoclonal antibody

Dilutions

1:1000

Preparation

Antigen: A synthetic phosphopeptide corresponding to residues surrounding Tyr516 of human TrkB

Formulation


Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Validations



Western blot detection of Phospho-TrkA (tyr490)/TrkB (Tyr516) in K562, C6, 3T3 cell lysates using Phospho-TrkA (tyr490)/TrkB (Tyr516) Rabbit mAb(1:1000 diluted). Predicted band size: 92kDa. Observed band size: 140kDa.