

Mouse IFNA5/IFNaG (His Tag) recombinant protein



Catalog Number: 502429

General Information

Gene Name Synonym

Interferon alpha-5; MCG123694

Protein Construction

A DNA sequence encoding the mouse IFNA5 (NP_034635.2) (Met1-Glu189) was expressed with a C-terminal polyhistidine tag.

Organism

Mouse

Expression Host

Human Cells

QC Testing

Purity

> 95 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per µg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Cys 24

Molecular Mass

The recombinant mouse IFNA5 comprises 177

amino acids and has a predicted molecular mass of 20.4 kDa. The apparent molecular mass of the protein is approximately 23.8 and 22.4 kDa in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Mouse IFNA5 / IFNaG Protein (His Tag) SDS-PAGE