

Anti-CXCR6 antibody



Catalog Number: 176429

Product name

Anti-CXCR6 antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to CXCR6

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from the n-terminal region of human cxcr6. aa range:1-50

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P

Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-1:300

Validations




Fig1:: Western Blot analysis of HeLa cells using Bonzo Polyclonal Antibody.. Secondary antibody [catalog#: HA1001) was diluted at 1:20000




Fig2:: Immunohistochemical analysis of paraffin-embedded human-pancreas, antibody was diluted at 1:100