

TSG-6 antibody



Catalog Number: 116429

Product name

TSG-6 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

TSG-6 Rabbit Polyclonal antibody. Positive WB detected in NIH/3T3 cells, A549 cells. Observed molecular weight by Western-blot: 39 kDa

Preparation

This antibody was obtained by immunization of TSG-6 recombinant protein (Accession Number: NM_007115). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



NIH/3T3 cells were subjected to SDS PAGE followed by western blot with Catalog No:116429(TNFAIP6 Antibody) at dilution of 1:600