

# MPO antibody



Catalog Number: 107429

## Product name

MPO antibody

## Specificity

Human; other species not tested.

## Antibody description

MPO Mouse Monoclonal antibody. Positive WB detected in HL-60 cells. Positive IHC detected in human tonsillitis tissue, human liver tissue. Observed molecular weight by Western-blot: 90 kDa

## Preparation

This antibody was obtained by immunization of MPO recombinant protein (Accession Number: NM\_000250). Purification method: Protein A purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgA

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

## Validations



HL-60 cells were subjected to SDS PAGE followed by western blot with Catalog No:107429(MPO Antibody) at dilution of 1:1000



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using (MPO Antibody) at dilution of 1:200