

Phospholipase C beta 1 antibody



Catalog Number: 113829

Product name

Phospholipase C beta 1 antibody

Immunogen

[Human PLCB1 Recombinant protein \(GST tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

Phospholipase C beta 1 Rabbit Polyclonal antibody. Positive WB detected in rat brain tissue. Observed molecular weight by Western-blot: 139 kDa

Preparation

This antibody was obtained by immunization of Phospholipase C beta 1 recombinant protein (Accession Number: NM_015192). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:500-1:5000

Validations



rat brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:113829(PLCB1 Antibody) at dilution of 1:1000