

# WIPI1 antibody



Catalog Number: 116829

## Product name

WIPI1 antibody

## Immunogen

[Human WIPI1 Recombinant protein \(GST tag\)](#)

## Specificity

Human, Mouse; other species not tested.

## Antibody description

WIPI1 Rabbit Polyclonal antibody. Positive WB detected in SW 1990 cells, HEK-293 cells, NIH/3T3 cells. Observed molecular weight by Western-blot: 49 kDa

## Preparation

This antibody was obtained by immunization of WIPI1 recombinant protein (Accession Number: NM\_017983). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:500-1:5000

## Validations



SW 1990 cells were subjected to SDS PAGE followed by western blot with Catalog No:116829(WIPI1 Antibody) at dilution of 1:1000