

# Anti-EPHX2 Rabbit antibody



Catalog Number: 167829

## Product name

Anti-EPHX2 Rabbit antibody

## Specificity

Human, Mouse

## Antibody description

EPHX2 Rabbit monoclonal antibody

## Preparation

Antigen: A synthetic peptide of human EPHX2

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,  
40%Glycerol, 0.01% sodium azide and 0.05% BSA

## Storage

Store at 4°C short term. Store at -20°C long term.  
Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications


WB, FC

## Dilutions

WB: 1/2000-1/10000

FC: 1/20

## Validations



Western blot detection of EPHX2 in A549,HL-60,U251,U87-MG cell lysates using EPHX2 Rabbit mAb(1:1000 diluted).Predicted band size:63kDa.Observed band size:63kDa.