

Dynamin 1 antibody

Catalog Number: 110129

Product name

Dynamin 1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

Dynamin 1 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, NIH/3T3 cells. Positive IP detected in mouse brain tissue. Observed molecular weight by Western-blot: 100 kDa

Preparation

This antibody was obtained by immunization of Dynamin 1 recombinant protein (Accession Number: XM_017014371). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IP: 1:1000-1:10000

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:110129(DNM1 antibody) at dilution of 1:2000



IP Result of anti-DNM1 (IP:Catalog No:110129, 4ug; Detection:Catalog No:110129 1:2000) with mouse brain tissue lysate 4000ug.