

# Neurokinin-1 receptor antibody



Catalog Number: 113129

## Product name

Neurokinin-1 receptor antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

Neurokinin-1 receptor Rabbit Polyclonal antibody. Positive WB detected in SH-SY5Y cells, mouse pancreas tissue. Observed molecular weight by Western-blot: 48-55 kDa

## Preparation

This antibody was obtained by immunization of Neurokinin-1 receptor recombinant protein (Accession Number: NM\_001058). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with Catalog No:113129(TACR1 Antibody) at dilution of 1:600