

# TIMM10 antibody



Catalog Number: 116129

## Product name

TIMM10 antibody

## Immunogen

[Human TIMM10 Recombinant protein \(GST tag\)](#)

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

TIMM10 Rabbit Polyclonal antibody. Positive IHC detected in human colon cancer tissue. Positive WB detected in human heart tissue, human liver tissue. Observed molecular weight by Western-blot: 10 kDa

## Preparation

This antibody was obtained by immunization of TIMM10 recombinant protein (Accession Number: NM\_012456). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

## Validations



human heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:116129(TIMM10 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:116129(TIMM10 antibody) at dilution of 1:50 (under 25x lens)