

Anti-GAB2 Mouse antibody



Catalog Number: 169129

Product name

Anti-GAB2 Mouse antibody

Specificity

Human

Antibody description

GAB2 Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant fragment of human GAB2 (AA: 332-476) expressed in E. Coli.

Formulation

Purified antibody in PBS with 0.05% sodium azide.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal


Ig Type

Mouse IgG2a


Applications

IHC-P, ICC/IF, FC


Validations



Western blot analysis using GAB2 mAb against HEK293 (1) and GAB2 (AA



Immunofluorescence analysis of HeLa cells using




GAB2 mouse mAb (green). Blue


Western blot analysis using GAB2 mAb against human GAB2 (AA

Anti-GAB2 Mouse antibody

Catalog Number: 169129



Immunohistochemical analysis of paraffin-embedded esophageal cancer tissues using GAB2 mouse mAb with DAB staining.



Flow cytometric analysis of HeLa cells using GAB2 mouse mAb (green) and negative control (red).

Immunohistochemical analysis of paraffin-embedded cervical cancer tissues using GAB2 mouse mAb with DAB staining.

