

HOXA6 antibody



Catalog Number: 111529

Product name

HOXA6 antibody

Specificity

Human; other species not tested.

Antibody description

HOXA6 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue.

Preparation

This antibody was obtained by immunization of HOXA6 recombinant protein (Accession Number: NM_024014). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:111529(HOXA6 Antibody) at dilution of 1:600