

General Information

Gene Name Synonym

GP1,2; GP1; GP2; GP2-delta

Protein Construction

A DNA sequence encoding the Zaire ebolavirus (strain Mayinga 1976) GP (AAC54887.1) (Met1-Gln650) was expressed with a polyhistidine tag at the C-terminus.

Organism

EBOV

Expression Host

Baculovirus-Insect Cells

QC Testing

Purity

> 95% as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg protein as determined by the LAL method.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Ile

Molecular Mass

The recombinant Zaire ebolavirus (strain Mayinga 1976) GP consists 629 amino acids and predicts a molecular mass of 69.3 kDa.

Formulation

Lyophilized from sterile 20mM Tris,500mM NaCl,10% glycerol, pH 8.0.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Ebola virus EBOV (subtype Zaire, strain Mayinga 1976) Glycoprotein / GP Protein (His Tag)