

SPAG11A antibody



Catalog Number: 115529

Product name

SPAG11A antibody

WB: 1:500-1:5000

IHC: 1:20-1:200


Specificity

Human; other species not tested.

Validations

Antibody description

SPAG11A Rabbit Polyclonal antibody. Positive WB detected in human testis tissue, human brain tissue. Positive IHC detected in human testis tissue, human kidney tissue, human lung tissue, human ovary tissue, human placenta tissue, human skin tissue, human spleen tissue. Observed molecular weight by Western-blot: 20kd




Preparation

This antibody was obtained by immunization of SPAG11A recombinant protein (Accession Number: BC058833). Purification method: Antigen affinity purified.

human testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:115529(SPAG11A antibody) at dilution of 1:300

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB, IHC

Dilutions

Recommended Dilution:

Immunohistochemical of paraffin-embedded human testis using Catalog No:115529(SPAG11A antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded

SPAG11A antibody



Catalog Number: 115529

human testis using Catalog No:115529(SPAG11A antibody) at dilution of 1:50 (under 40x lens)