

Phospho-BCAR1-Y165 pAb



Catalog Number: 166529

Product name

Phospho-BCAR1-Y165 pAb

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to Phospho-BCAR1-Y165

Preparation

Antigen: A phospho specific peptide corresponding to residues surrounding Y165 of human BCAR1

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB 1:500 - 1:2000

Validations



Western blot - Phospho-BCAR1-Y165 pAb

Western blot analysis of extracts of PC-12 cells, using Phospho-BCAR1-Y165 antibody at 1:1000 dilution. PC-12 cell lysate were treated by CIP (20ul CIP for each 400ul cell lysate) at 37u2103 for 1 hour .Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution.Lysates/proteins: 25ug per lane.Blocking buffer: 3% BSA.Detection: ECL Basic Kit .Exposure time: 90s.