

RRP7A antibody



Catalog Number: 114929

Product name

RRP7A antibody

IHC: 1:20-1:200

IF: 1:10-1:100

Specificity

Human, Mouse; other species not tested.

Antibody description

RRP7A Rabbit Polyclonal antibody. Positive IF detected in A431 cells. Positive IHC detected in human brain tissue, human stomach cancer tissue. Positive WB detected in mouse thymus tissue. Observed molecular weight by Western-blot: 42 kDa

Preparation

This antibody was obtained by immunization of RRP7A recombinant protein (Accession Number: NM_015703). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IF, IHC

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


Validations



mouse thymus tissue were subjected to SDS PAGE followed by western blot with Catalog No:114929(RRP7A Antibody) at dilution of 1:1000



Immunofluorescent analysis of A431 cells using Catalog No:114929(RRP7A Antibody) at dilution of 1:25 and Rhodamine-Goat anti-Rabbit IgG



Immunohistochemistry of paraffin-embedded human brain tissue slide using Catalog

RRP7A antibody



Catalog Number: 114929

No:114929(RRP7A Antibody) at dilution of 1:50 (under 40x lens)