

Anti-phospho-EEF2/Elongation factor 2 (Thr56/Thr58) Rabbit antibody



Catalog Number: 167929

Product name

Monoclonal

Anti-phospho-EEF2/Elongation factor 2 (Thr56/Thr58) Rabbit antibody

Ig Type

Rabbit IgG

Specificity

Human, Mouse, Rat

Applications

WB

Antibody description

phospho-EEF2/Elongation factor 2 (Thr56/Thr58)
Rabbit monoclonal antibody


Dilutions

WB: 1/1000

Preparation

Antigen: A synthetic phosphopeptide corresponding to residues surrounding Thr56/Thr58 of human EEF2/Elongation factor 2

Validations



Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,
40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Western blot detection of EEF2/Elongation factor 2 (phospho-Thr56/Thr58) in HeLa,HL-60 using EEF2/Elongation factor 2 (phospho-Thr56/Thr58) Rabbit mAb(1:1000 diluted)

Clonality