

# Anti-hnRNP K Rabbit antibody



Catalog Number: 168929

## Product name

Anti-hnRNP K Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

hnRNP K Rabbit monoclonal antibody

## Preparation

Antigen: A synthetic peptide of human hnRNP K

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications


WB, IHC-P

## Dilutions

1:5000-1:10000

1:20-1:100

## Validations



Western blot detection of hnRNP K in K562, Rat Brain, C6, 3T3, Hela cell lysates using hnRNP K Rabbit mAb(1:1000 diluted). Predicted band size: 51kDa. Observed band size: 62kDa.