

DHX29 antibody



Catalog Number: 109929

Product name

DHX29 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

DHX29 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, HepG2 cells, HT-1080 cells, K-562 cells, MCF7 cells, PC-3 cells.

Observed molecular weight by Western-blot:
150-155 kDa

Preparation

This antibody was obtained by immunization of DHX29 recombinant protein (Accession Number: NM_019030). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol
pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:109929(DHX29 antibody) at dilution of 1:500