

Anti-OR7C2 antibody



Catalog Number: 176229

Product name

Anti-OR7C2 antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to OR7C2

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human or7c2. aa range:198-247

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB

Dilutions

WB: 1:500-1:2000

Validations



Fig1.; Western Blot analysis of various cells using Olfactory receptor 7C2 Polyclonal Antibody