

Anti-WDR77 (8A10) Mouse antibody



Catalog Number: 168229

Product name

Anti-WDR77 (8A10) Mouse antibody

Specificity

Human, Mouse, Rat

Antibody description

WDR77 (8A10) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant human WDR77 protein fragments expressed in E.coli.

Formulation

0.02% NaN₃ 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications


WB, ICC/IF

Dilutions


WB: 1/1000

;ICC: 1/100

Validations



Immunocytochemistry staining of HeLa cells fixed in 1% Paraformaldehyde and then permeabilized in 0.1% Triton X-100, next using WDR77 mouse



mAb (dilution 1:100).

Western blot detection of WDR77 in C6, 3T3 and K562 cell lysates using WDR77 mouse mAb (1:1000 diluted). Predicted band size: 42KDa. Observed band size: 42KDa.