

Anti-CHI3L1 antibody



Catalog Number: 176629

Product name

Anti-CHI3L1 antibody

Specificity

Human

Antibody description

Mouse monoclonal antibody to CHI3L1

Preparation

This antigen of this antibody was recombinant protein within extracellular domain of human chi3l1.

Formulation

Liquid, 1*TBS (pH7.4), 0.2% BSA, 50% Glycerol.
Preservative: 0.05% Sodium Azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C. Avoid repeated freeze / thaw cycles.

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications

WB

Dilutions

WB:1:500-1:2,000

Validations




Fig1.; Western blot analysis of YKL-40 on THP-1 cell lysates. Proteins were transferred to a PVDF membrane and blocked with 5% NBTM/TBST for 1 hour at room temperature. The primary antibody (1/500) was used in 5% NBTM/TBST at room temperature for 2 hours. Goat Anti-Mouse IgG - HRP Secondary Antibody (HA1006) at 1:5,000 dilution was used for 1 hour at room temperature.