

Human ALDOA/ALDA (His tag) recombinant protein



Catalog Number: 519629

General Information

Gene Name Synonym

Fructose-Bisphosphate Aldolase A, Lung Cancer Antigen NY-LU-1, Muscle-Type Aldolase, ALDOA, ALDA

Protein Construction

A DNA sequence encoding the Human ALDOA/ALDA Pro2-Tyr364 is expressed with a 6His tag at the C-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μ g protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 μ m filtered solution of 20mM TrisHCl, 100mM NaCl, 20% Glycerol, pH 8.0.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Human ALDOA/ALDA (His tag) recombinant protein