

Entactin antibody

Catalog Number: 110329

Product name

Entactin antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

Entactin Rabbit Polyclonal antibody. Positive WB detected in human placenta tissue, mouse heart tissue. Positive IP detected in human placenta tissue. Observed molecular weight by Western-blot: 136 kDa

Preparation

This antibody was obtained by immunization of Entactin recombinant protein (Accession Number: BC045606). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:1000

Validations



human placenta tissue were subjected to SDS PAGE followed by western blot with Catalog No:110329(NID1 antibody) at dilution of 1:200



IP Result of anti-NID1 (IP:Catalog No:110329, 4ug; Detection:Catalog No:110329 1:200) with human placenta tissue lysate 4000ug.