

# TRAM2 antibody



Catalog Number: 116329

## Product name

TRAM2 antibody

## Immunogen

[Human TRAM2 Recombinant protein \(GST tag\)](#)

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

TRAM2 Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells, Sp2/0 cells. Observed molecular weight by Western-blot: 43kd,66kd

## Preparation

This antibody was obtained by immunization of TRAM2 recombinant protein (Accession Number: NM\_012288). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:200-1:2000

## Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:116329(TRAM2 antibody) at dilution of 1:500