

# C4BPB antibody



Catalog Number: 108729

## Product name

Polyclonal

C4BPB antibody

## Ig Type

## Specificity

Rabbit IgG

Human, Mouse; other species not tested.

## Applications

## Antibody description

ELISA, WB

C4BPB Rabbit Polyclonal antibody. Positive WB detected in human plasma tissue, HeLa cells, mouse ovary tissue. Observed molecular weight by Western-blot: 45 kDa

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Preparation

This antibody was obtained by immunization of C4BPB recombinant protein (Accession Number: NM\_001017364). Purification method: Antigen affinity purified.

## Validations


## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality



human plasma tissue were subjected to SDS PAGE followed by western blot with Catalog No:108729(C4BPB Antibody) at dilution of 1:600