

CTPS2 antibody

Catalog Number: 109729

Product name

CTPS2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

CTPS2 Rabbit Polyclonal antibody. Positive IHC detected in human prostate cancer tissue, human colon tissue. Positive WB detected in mouse testis tissue, human adrenal gland tissue. Observed molecular weight by Western-blot: 66 kDa

Preparation

This antibody was obtained by immunization of CTPS2 recombinant protein (Accession Number: NM_019857). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:109729(CTPS2 antibody) at dilution of 1:500



Immunohistochemistry of paraffin-embedded human prostate cancer tissue slide using Catalog No:109729(CTPS2 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human prostate cancer tissue slide using Catalog No:109729(CTPS2 Antibody) at dilution of 1:50

CTPS2 antibody



Catalog Number: 109729

(under 40x lens)