

DTX3L antibody



Catalog Number: 110039

Product name

DTX3L antibody

Specificity

human; other species not tested.

Antibody description

DTX3L Rabbit Polyclonal antibody. Positive WB detected in human lung tissue. Observed molecular weight by Western-blot: 84kd

Preparation

This antibody was obtained by immunization of DTX3L recombinant protein (Accession Number: NM_138287). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



human lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:110039(DTX3L antibody) at dilution of 1:500