

GLTP antibody



Catalog Number: 111039

Product name

GLTP antibody

Immunogen

[Human GLTP Recombinant protein \(GST tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

GLTP Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, A2780 cells. Observed molecular weight by Western-blot: 24 kDa

Preparation

This antibody was obtained by immunization of GLTP recombinant protein (Accession Number: NM_016433). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:200-1:2000

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:111039(GLTP antibody) at dilution of 1:500