

KHK antibody

Catalog Number: 112039

Product name

KHK antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

KHK Rabbit Polyclonal antibody. Positive WB detected in human liver tissue, mouse liver tissue. Positive IP detected in mouse liver tissue. Observed molecular weight by Western-blot: 33-35 kDa

Preparation

This antibody was obtained by immunization of KHK recombinant protein (Accession Number: NM_000221). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

Validations



human liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:112039(KHK antibody) at dilution of 1:1200



IP Result of anti-KHK (IP:Catalog No:112039, 4ug; Detection:Catalog No:112039 1:500) with mouse liver tissue lysate 4000ug.