

Anti-FGFR3 Rabbit antibody



Catalog Number: 168039

Product name

Anti-FGFR3 Rabbit antibody

Specificity

Human

Antibody description

FGFR3 Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human FGFR3

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, FC

Dilutions


WB: 1/2000

IHC: 1/50

ICC/IF: 1/100

FC: 1/20

Validations



Western blot detection of FGFR3 in HeLa cell lysates using FGFR3 Rabbit mAb(1:1000 diluted). Predicted band size:88kDa.Observed band size:88kDa.