

Anti-HMG4 Rabbit antibody



Catalog Number: 169039

Product name

Anti-HMG4 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

HMG4 Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human HMG4

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

1:2000-1:10000

Validations



Western blot detection of HMG4 in K562,C6,3T3,HeLa cell lysates using HMG4 Rabbit mAb(1:1000 diluted).Predicted band size:23kDa.Observed band size:23kDa.