

# Human DFFA/DFF1/DFF45/H13 (His tag) recombinant protein



Catalog Number: 519039

## General Information

### Gene Name Synonym

DNA Fragmentation Factor Subunit Alpha, DNA Fragmentation Factor 45 kDa Subunit, DFF-45, Inhibitor of CAD, ICAD, DFFA, DFF1, DFF45, H13

### Protein Construction

A DNA sequence encoding the Human DFFA/DFF1/DFF45/H13 Met1-Thr331 is expressed with a 6His tag at the N-terminus.

### Organism

Human

### Expression Host

E. coli

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per µg protein as determined by the LAL

method.

### Stability

Store for up to 6 months at -20°C

### Formulation


Supplied as a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH7.4.

## Usage Guide

### Storage

Store at <-20°C, minimize freeze-thaw cycles.

### SDS-PAGE



Human DFFA/DFF1/DFF45/H13 (His tag)  
recombinant protein