

# ATRIP Polyclonal Antibody



Catalog Number: 164439

## Product name

ATRIP Polyclonal Antibody

## Specificity

Human

## Antibody description

Polyclonal antibody to ATRIP

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human ATRIP (NP\_569055.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IF, IP


## Dilutions

WB 1:500 - 1:2000

IF 1:50 - 1:200

IP 1:20 - 1:50

## Validations



## Immunoprecipitation - ATRIP Polyclonal Antibody

Immunoprecipitation analysis of 200ug extracts of HeLa cells using 1ug ATRIP antibody . Western blot was performed from the immunoprecipitate using ATRIP antibody at a dilution of 1:1000.