

OTUD7B antibody



Catalog Number: 107439

Product name

OTUD7B antibody

Specificity

Human; other species not tested.

Antibody description

OTUD7B Mouse Monoclonal antibody. Positive WB detected in Jurkat cells. Observed molecular weight by Western-blot: 92 kDa

Preparation

This antibody was obtained by immunization of OTUD7B recombinant protein (Accession Number: BC020622). Purification method: Protein A purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG2b

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:107439(OTUD7B Antibody) at dilution of 1:1000