

Anti-TFIIB (2F6) Mouse antibody



Catalog Number: 168439

Product name

Anti-TFIIB (2F6) Mouse antibody

Ig Type

Mouse IgG2b

Specificity

Human, Mouse, Monkey

Applications

WB

Antibody description

TFIIB (2F6) Mouse monoclonal antibody

Dilutions

WB: 1/500


Preparation

Antigen: Purified recombinant human TFIIB protein fragments expressed in E.coli.

Validations

Formulation

Ascites containing with 0.03% Proclin300 and 50% glycerol.



Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Western blot detection of TFIIB in Jurkat, MCF7, COS7, 3T3 and HeLa cell lysates using TFIIB mouse mAb (1:500 diluted). Observed band size: 34 kDa

Clonality

Monoclonal