

Collagen Type VI antibody



Catalog Number: 109439

Product name

Collagen Type VI antibody

IHC: 1:20-1:200


Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

Collagen Type VI Rabbit Polyclonal antibody.
Positive WB detected in human heart tissue.
Positive IHC detected in human heart tissue.
Observed molecular weight by Western-blot: 140 kDa

Validations




Preparation

This antibody was obtained by immunization of Collagen Type VI recombinant protein (Accession Number: XM_011529451). Purification method: Antigen affinity purified.

human heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:109439(COL6A2 antibody) at dilution of 1:500

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Immunohistochemical of paraffin-embedded human heart using Catalog No:109439(COL6A2 antibody) at dilution of 1:50 (under 40x lens)

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB, IHC

Dilutions

Recommended Dilution:

WB: 1:200-1:2000



Immunohistochemical of paraffin-embedded human heart using Catalog No:109439(COL6A2 antibody) at dilution of 1:50 (under 10x lens)