

Mouse Resistin / RETN (Fc Tag) recombinant protein



Catalog Number: 502839

General Information

Gene Name Synonym

Adipose tissue-specific secretory factor; Adipose-specific cysteine-rich secreted protein A12-alpha; Cysteine-rich secreted protein FIZZ3

Protein Construction

A DNA sequence encoding the mouse RETN (Q99P87) (Ser 21-Ser 114) was fused with the Fc region of human IgG1 at the N-terminus.

Organism

Mouse

Expression Host

Human Cells

QC Testing

Purity

> 80 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Glu

Molecular Mass

The recombinant mouse RETN/Fc is a disulfide-

linked homodimer. The reduced monomer comprises 354 amino acids and has a calculated molecular mass of 38.6 kDa. The apparent molecular mass of the monomer is approximately 40 kDa in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide

Storage


Store it under sterile conditions at -20°C to -80°C .

It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Mouse Resistin / ADFP / RETN Protein (Fc Tag)
SDS-PAGE