

# SUN1 antibody

Catalog Number: 115839

## Product name

SUN1 antibody

## Specificity

Human; other species not tested.

## Antibody description

SUN1 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, A431 cells, Jurkat cells, K-562 cells. Positive IP detected in Jurkat cells. Observed molecular weight by Western-blot: 90-100 kDa

## Preparation

This antibody was obtained by immunization of SUN1 recombinant protein (Accession Number: NM\_001171946). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:500-1:5000

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:115839(UNC84A Antibody) at dilution of 1:1000



IP Result of anti-UNC84A (IP:Catalog No:115839, 4ug; Detection:Catalog No:115839 1:1000) with Jurkat cells lysate 3200ug.