

Product name

AGTRAP antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

AGTRAP Rabbit Polyclonal antibody. Positive IF detected in HEK-293 cells. Positive IHC detected in human pancreas cancer tissue.

Preparation

This antibody was obtained by immunization of AGTRAP recombinant protein (Accession Number: NM_001040194). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, IF


Dilutions

Recommended Dilution:


IHC: 1:10-1:100

IF: 1:20-1:200

Validations



Immunohistochemical of paraffin-embedded human pancreas cancer using Catalog No:107839(AGTRAP antibody) at dilution of 1:10 (under 10x lens)



Immunofluorescent analysis of (-20oc Ethanol) fixed HEK-293 cells using Catalog No:107839(AGTRAP Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)