

DDX56 antibody



Catalog Number: 109839

Product name

DDX56 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

DDX56 Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells, mouse brain tissue. Observed molecular weight by Western-blot: 65 kDa

Preparation

This antibody was obtained by immunization of DDX56 recombinant protein (Accession Number: NM_019082). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:109839(DDX56 Antibody) at dilution of 1:1000