

PPM2C antibody

Catalog Number: 114139

Product name

PPM2C antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

PPM2C Rabbit Polyclonal antibody. Positive WB detected in mouse skeletal muscle tissue, mouse brain tissue. Positive IP detected in mouse brain tissue. Observed molecular weight by Western-blot: 53 kDa

Preparation

This antibody was obtained by immunization of PPM2C recombinant protein (Accession Number: XM_047421906). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IP: 1:200-1:1000

Validations



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:114139(PPM2C Antibody) at dilution of 1:300



IP Result of anti-PPM2C (IP:Catalog No:114139, 4ug; Detection:Catalog No:114139 1:300) with mouse brain tissue lysate 5200ug.