

Anti-OR10A7 antibody



Catalog Number: 176139

Product name

Anti-OR10A7 antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to OR10A7

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human or10a7. aa range:234-283

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB

Dilutions

WB: 1:500-1:2000

Validations




Fig1.; Western Blot analysis of K562 cells using Olfactory receptor 10A7 Polyclonal Antibody