

# Anti-Cytokeratin 18 (2C3) Mouse antibody



Catalog Number: 168139

## Product name

Anti-Cytokeratin 18 (2C3) Mouse antibody

## Dilutions

WB: 1/500-1000

## Specificity

Human, Mouse

## Antibody description

Cytokeratin 18 (2C3) Mouse monoclonal antibody

## Preparation

Antigen: Recombinant Protein of Keratin type I cytoskeletal 18

## Formulation

PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal


## Ig Type

Mouse IgG1


## Applications

WB

## Validations



Western blot detection of Keratin in Goat Muscle, CHO-K1, C6, HeLa, 3T3 cell lysates using Keratin 18 (2C3) Mouse mAb (1:1000 diluted). Predicted band size: 46 kDa. Observed band size: 46 kDa.



Western blot analysis of 1) HeLa, 2) HepG2, 3) Mouse Skeletal Muscle tissue, (2F7) diluted at 1:1000.