

# MGAT4C antibody



Catalog Number: 112539

## Product name

MGAT4C antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

MGAT4C Rabbit Polyclonal antibody. Positive WB detected in SH-SY5Y cells, mouse testis tissue. Positive IP detected in mouse testis tissue. Observed molecular weight by Western-blot: 55 kDa

## Preparation

This antibody was obtained by immunization of MGAT4C recombinant protein (Accession Number: XM\_017019150). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:1000

## Validations



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with Catalog No:112539(MGAT4C antibody) at dilution of 1:300



IP Result of anti-MGAT4C (IP:Catalog No:112539, 4ug; Detection:Catalog No:112539 1:400) with mouse testis tissue lysate 4000ug.