

RANBP6 antibody



Catalog Number: 114539

Product name

RANBP6 antibody

Specificity

human, mouse; other species not tested.

Antibody description

RANBP6 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, mouse skeletal muscle tissue. Observed molecular weight by Western-blot: 125kd

Preparation

This antibody was obtained by immunization of RANBP6 recombinant protein (Accession Number: NM_012416). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:114539(RANBP6 antibody) at dilution of 1:200