

# UBE2V1 antibody



Catalog Number: 116539

## Product name

UBE2V1 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

UBE2V1 Rabbit Polyclonal antibody. Positive WB detected in human spleen tissue, BxPC-3 cells, mouse spleen tissue, Transfected HEK-293 cells. Observed molecular weight by Western-blot: 20 kDa

## Preparation

This antibody was obtained by immunization of UBE2V1 recombinant protein (Accession Number: NM\_001032288). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:1000

## Validations



human spleen tissue were subjected to SDS PAGE followed by western blot with Catalog No:116539(UBE2V1 antibody) at dilution of 1:300