

Mouse Mmp12/Mme/Mmel (His tag) recombinant protein



Catalog Number: 518539

General Information

Gene Name Synonym

Macrophage metalloelastase, MME, Matrix metalloproteinase-12, MMP-12, Mmp12, Mme, Mmel

Protein Construction

A DNA sequence encoding the Mouse Mmp12/Mme/Mmel Ala29-Cys473 is expressed with a 6His tag at the C-terminus.

Organism

Mouse

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of PBS, pH 7.4.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Mouse Mmp12/Mme/Mmel (His tag) recombinant protein